



1 Harrier Court  
 Airfield Business Park  
 Elvington  
 York YO41 4AU  
 United Kingdom  
 Tel: +44 (0) 1904 607600  
 Fax: +44 (0) 1904 607601  
 Web: www.vet-way.com

## BIOSHIELD IIIR

Product	Code	Case Quantity	Ext Dim. L x W x H (mm)	Volumetric Weight (kg)	Total Weight (kg)
BIOSHIELD IIIR Box	BSE043	1	470 x 500 x 340	13.31	0.85
BIOSHIELD IIIR Carcass Kit	PPE001	1	N/A	N/A	0.13

### Components Supplied

- Outer box printed in compliance with transport IATA (air) and ADR (road) regulation
- Pathoseal Liquitite 20/32 leak-proof bags
- Absorbent material

### Features

- UN certified for the shipment of diagnostic specimens
- Flat-packed for easy storage
- Easy to assemble
- Manufactured from strong, durable and high quality Correx plastic
- Complete with handles to encourage careful and easy handling

The Bioshield IIIR carcass shipper is to be used for the transportation of solid, diagnostic (UN3373), samples.

The Bioshield IIIR Box needs to be used in conjunction with the Bioshield IIIR Carcass Kit for this system to be compliant. Each Bioshield IIIR Carcass Kit consists of 2 x 95kPa leak-proof bags and absorbent. Samples are to be placed inside both 95kPa leak-proof bags with the absorbent provided and then placed inside the Bioshield IIIR outer Correx box for transportation.

Should multiple samples need to be shipped together in the same Bioshield IIIR outer, a Bioshield IIIR Carcass Kit will need to be purchased for each sample being shipped.

