



1 Harrier Court
 Airfield Business Park
 Elvington
 York YO41 4AU
 United Kingdom
 Tel: +44 (0) 1904 607600
 Fax: +44 (0) 1904 607601
 Web: www.vet-way.com

BIOTHERM 5 & 10

| Product | Code | Case Quantity | Volumetric Wt. (kg) | +2°C to +8°C Validated | Average Temp (°C) During Test | Dry Ice (-20°C) Performance Hours |
|-----------------------|--------------------------|-------------------|-------------------------|--------------------------|-------------------------------|-----------------------------------|
| Biotherm 5 | BT002 | 1 | 3.4 | TBC | 20 | 72 |
| Dry Ice Capacity (kg) | Dry Ice Loss Rate (gm/h) | Dry Ice Used (kg) | Ext Dim. L x W x H (mm) | Int. Dim. L x W x H (mm) | Internal Vol. (Ltr) | Total Wt. (kg) |
| 6.8 | 45 | 3.5 | 270 x 270 x 275 | 170 x 170 x 170 | 4.77 | 0.81 |

| Product | Code | Case Quantity | Volumetric Wt. (kg) | +2°C to +8°C Validated | Average Temp (°C) During Test | Dry Ice (-20°C) Performance Hours |
|-----------------------|--------------------------|-------------------|-------------------------|--------------------------|-------------------------------|-----------------------------------|
| Biotherm 10 | BT001 | 1 | 6.05 | 72 | 20 | 96 |
| Dry Ice Capacity (kg) | Dry Ice Loss Rate (gm/h) | Dry Ice Used (kg) | Ext Dim. L x W x H (mm) | Int. Dim. L x W x H (mm) | Internal Vol. (Ltr) | Total Wt. (kg) |
| 8.5 | 70 | 6.8 | 360 x 315 x 320 | 210 x 250 x 185 | 9.71 | 1.18 |





1 Harrier Court
Airfield Business Park
Elvington
York YO41 4AU
United Kingdom
Tel: +44 (0) 1904 607600
Fax: +44 (0) 1904 607601
Web: www.vet-way.com

BIO THERM 5 & 10

Components Supplied

- BIO THERM outer box printed in compliance with transport regulations IATA (air) and ADR (road)
- Plastic handle
- EPS insulated shipper
- Infectious specimen label
- Dry ice label
- Overpack label
- Security seal label

Features

- Re-usable, low cost and highly effective
- Ideal for shipping a wide range of dry ice or temperature controlled specimens
- Carry handle encourages careful transportation
- Temperature tested performance reports available on request
- Validated -2°C to -8°C Coolpac kits available
- Supplied complete ready to ship
- Can be used with DGP PATHOPAK products for diagnostic and infectious shipments

The BIO THERM 5 & 10 are insulated shippers to transport samples at a controlled temperature.

These insulated shippers are ideal to be used in combination with the PATHOPAK 1Ltr (BIO THERM 10) and the PATHOPOUCH A4 (BIO THERM 5), to keep blood or tissue between 2°C to 8°C in temperature.

The BIO THERM 5 & 10 require the Summer or Winter Temperature Control Validated COOLPAC Kits to be purchase separately.

Specimens are placed in the ambient specimen shippers (PATHOPAK 1Ltr or PATHOPOUCH A4) and then placed inside the applicable insulated shipper with the validated COOLPAC kits.

