



1 Harrier Court
Airfield Business Park
Elvington
York YO41 4AU
United Kingdom
Tel: +44 (0) 1904 607600
Fax: +44 (0) 1904 607601
Web: www.vet-way.com



PATHOPAK Group

Product	Code	Case Quantity	Ext Dim. L x W x H (mm)	Int. Dim. Diam x H (mm)	Opening Diam (mm)	Volumetric Wt (kg)	Total Wt. (kg)
PATHOPAK 1Ltr	PP001	20	105 x 105 x 185	80 x 155	70	0.34	0.20
DGP Bio-Bottle 2Ltr	BB002	12	125 x 125 x 245	110 x 205	87	0.64	0.38
DGP Bio-Bottle 3Ltr	BB003	12	150 x 150 x 240	105 x 213	120	0.9	0.41

Components Supplied

- 1L leak-proof polypropylene container with 'o' ring lid (95kPa IATA test compliant)
- Outer box printed in compliance with transport regulations IATA (air) and ADR (road)
- Absorbent material with 50ml capacity
- Bubble wrap pouch
- Infectious class label
- Security Seal label

Features

- Supplied complete ready to ship
- UN certified for infectious/ diagnostic specimens
- Strong, durable and highly
- Protective
- Re-usable
- Cost-effective
- Easy to use
- Ideal for shipping a wide range of primary containers including blood tubes, vials swabs universal ports and petri dishes

The PATHOPAK 1Ltr, 2Ltr & 3Ltr are specimen shippers which can be used to send blood samples and swabs to the laboratory for both infectious (UN2814/2900) and diagnostic (UN3373) samples.

The PATHOPAK Group are ideal for shipping a wide variety of primary containers. A primary container must be used in conjunction with the components provided with these specimentshippers. The primary container with the specimen must be placed inside the container with the absorbent and bubble wrap bag provided. The 1Ltr, 2Ltr & 3Ltr container must then be placed inside the rigid and durable outer box.

