



1 Harrier Court  
Airfield Business Park  
Elvington  
York YO41 4AU  
United Kingdom  
Tel: +44 (0) 1904 607600  
Fax: +44 (0) 1904 607601  
Web: www.vet-way.com



## +2°C to +8°C Temperature Control Validated COOLPAC Kits

| Product              | Code  | Case Quantity | COOLPACS | EPS Spacers | Validated Hours | Total Wt. (kg) |
|----------------------|-------|---------------|----------|-------------|-----------------|----------------|
| BIOTHERM 5 (Summer)  | BT034 | 1             | 3        | 3           | 46              | 1.4            |
| BIOTHERM 10 (Summer) | BT013 | 1             | 3        | 3           | 76              | 1.6            |
| BIOTHERM 10 (Winter) | BT014 | 1             | 1        | 1           | 96              | 0.99           |

### Components Supplied

- As indicated by the table below

### Features

- +2°C to +8°C temperature control validated COOLPAC kits
- Temperature control devices should be pre-conditioned as instructed
- Should be chilled or frozen 24hrs prior to use
- Cost-effective
- To be used with BIOTHERM range (not supplied)

The Temperature Control Validated COOLPAC Kits are used with the BIOTHERM range for keeping samples between +2°C to +8°C in temperature during transportation.

The Validated COOLPAC Kits should be placed inside the insulated shippers and around the sample being contained in an ambient shipper.

