



1 Harrier Court
Airfield Business Park
Elvington
York YO41 4AU
United Kingdom
Tel: +44 (0) 1904 607600
Fax: +44 (0) 1904 607601
Web: www.vet-way.com

PATHOPAK 20Ltr

Product	Code	Case Quantity	Ext Dim. L x W x H (mm)	Int. Dim. Diam x H (mm)	Opening Diam (mm)	Volumetric Wt (kg)	Total Wt. (kg)
PATHOPAK 20Ltr	PP010	1	360 x 365 x 375	290 x 300	312	8.21	2.20

Components Supplied

- 20L leak-proof container (95kPa IATA test compliant)
- Heavy duty Correx outer box printed in compliance with transport regulations IATA (air) and ADR (road)
- Absorbent material
- Leak-proof liner
- Security Seal label

Features

- System tested to conform with PI650
- Easy access to keg reduces filling time
- Strong, durable and highly protective
- Re-usable
- Durable
- Cost-effective
- Ideal for shipping a wide range of primary containers including blood tubes, vials swabs universal ports and petri dishes

The PATHOPAK 20Ltr shipper can be used for the transportation of diagnostic (UN3373) specimens by road and air and is most suitable for transporting multiple liquid specimens.

To comply with regulations during transportation, specimens need to be placed within the leak-proof bag with the absorbent provided, within the 20Ltr container and then placed inside the outer Correx box. Alternatively, the PATHOPAK 20Ltr is ideal for shipping a wide range of primary containers, which can be placed inside the 20Ltr container, using the leak-proof liner as additional protection.

