



1 Harrier Court
Airfield Business Park
Elvington
York YO41 4AU
United Kingdom
Tel: +44 (0) 1904 607600
Fax: +44 (0) 1904 607601
Web: www.vet-way.com

PATHOPOUCH 95kPa A4 Pouch

Product	Code	Case Quantity	95kPa Test Compliant	Components Included	Opening Diam. (mm)	Size
PATHOPOUCH 95kPa (A4) c/w 6 Bay Absorbent	BP001	50	IATA (Air)	Absorbent	185	340 x 240

Components Supplied

- As indicated by the table below

Features

- 95kPa IATA test compliant for air and road transportation
- Easy to use tamper evident seal
- A4 size sequentially numbered
- Leak-proof
- Third party tested

The PATHOPOUCH 95kPa A4 pouch is used to transport blood, tissue samples and swabs to the laboratory.

This product must be used in conjunction with a rigid outer box and primary container and is ideal for being used in combination with the BIOTHERM 5 to send temperature controlled samples for analysis.

The sample is to be placed in a primary container within the pouch with the absorbent provided. The pouch is then placed within the BIOTHERM 5 temperature controlled specimen shipper.

